FOUR HOOKS HANGER FOR HYPODERMOCLYSIS PNT 698 IR



Dimensional Data 147 / 247 Dimensions (L x W x H) Cm Capacity max. Kg 28 Net weight Kg 20,5 Gross weight Kg 22,0 Volume in cube Mt. 0,06 Number of items 2 Packing (L x W x H) Cm (Base) 60 x 14 x 62 Packing (L x W x H) Cm (Rod) 10 x 10 x 150

Four hooks hanger for hypoderomoclysis, stainless steel rod for the hooksm adjustable height by a clamp.

Mounted on a 5 spoke aluminium base with swivel tyre castors D.50 mm (3 atistatic + 2 with brake)

Price without transport and taxes: 115,00 €

Contact us for a quotation!

Technical Data:

Structure:

Rod: stainless steel tubular frame AISI 304 ENG 10217/7 diameter 25 mm thickness 1,5, diameter 20 mm tickness 1,5 and diameter 18mm thickness 1. Base: die- cast aluminium.

Finish:

Rod: polished stainless steel Base: polished aluminium

Other components:

Such as screws and other mechanical parts, different materials suitable for the purpose are used.

Cleaning:

STAINLESS STEEL

Normal detergents can be used for cleaning, provided they do not contain chlorine and its compouds. When cleaning, it is important to follow the satin finish of stainless steel, without exerting too much pressure, For extremly stubborn stains, specific products for stainless steel should be used.

Important notes:

- 1. Do not use scouring pads or anything else that may scratch or scrape.
- 2. Avoid using muratic or bleach,
- 3. Do not leave iron objects on the tray for a long time.

ALUMINIUM

Clean with a slightly damp cloth and normal detergents, do not se sollvents or corrosive chemicals.

Telefon: +49 341 2612498 Online: www.cs-humanmobility.com